

Made in
the USA

SPANDAGRIP™

TUBULAR ELASTIC SUPPORT BANDAGE

SINGLE UNITS

LATEX FREE

SpandaGrip™ stays in place without pins, clips or tape. A convenient alternative to flat wrap bandages. Easy to apply. Not bulky, which allows freedom of movement and comfort.

Delivers consistent support with a single layer.



Order #	Description and Color	Size	Suggested Application
SAG25036	B - Natural	2-1/2" x 36"	Small hands and arms
SAG23436	C - Natural	2-3/4" x 36"	Medium arms, small ankles
SAG10336	D - Natural	3" x 36"	Large arms, med. ankles, small knees
SAG35036	E - Natural	3-1/2" x 36"	Large ankles, med. knees, small thighs
SAG10436	F - Natural	4" x 36"	Large knees, medium thighs
SAG45036	G - Natural	4-1/2" x 36"	Large thighs

PRODUCT INFORMATION

Perfect for the effective management of:

- Sprains and Strains
- Support following Cast Removal
- Soft-tissue Injuries
- Joints and Muscles
- General Edema
- Post-graft Pressure for Burn Patients

36 Inches Per Unit / 12 Units Per Case

M
MEDI-TECH
International Corporation

(800) 333-0109

www.medi-techintl.com