



# **SPANDAGRIP™**

## **TUBULAR ELASTIC SUPPORT BANDAGE**

### **LAUNDERING INSTRUCTIONS**

- **CAN BE WASHED UP TO 20 TIMES OR MORE WITHOUT THE LOSS OF ELASTICITY**
- **HAND WASH - WARM WATER WITH MILD LAUNDRY DETERGENT OR DISHWASHING LIQUID**
- **DRIP/HANG TO DRY**
- **DO NOT USE BLEACH**

Medi-Tech International Corp. Brooklyn, NY  
(800) 333-0109  
[www.medi-techintl.com](http://www.medi-techintl.com)