

Medi-Elbow & Heel Safeguard™

Skin Pressure Protection

Not Made with Natural Rubber Latex



Medi-Elbow & Heel Safeguard...
with foam padding decreases friction and shear on the elbows and heels.

Helps avoid dermal ulcers along with other dermal ulcer therapy by padding & shielding sensitive skin.

PRODUCT INFORMATION

Order #	Size / Color	Circumference & Style	Unit of Measure
MTICSM59	Small / Yellow	Up to 16" Regular	12 Pair / Case
MTICMED60	Medium / Blue	Up to 17" Regular	12 Pair / Case
MTICLG61	Large / White	Up to 19" Regular	12 Pair / Case
MTICXLG62	X-Large / Green	Up to 23" Regular	12 Pair / Case
MTICXXLG63	XX-Large / Beige	Up to 25" Regular	12 Pair / Case

- The soft, durable fabric consists of: 72% Acrylic - 20% Spandex - 8% Nylon.

- A comfortable layer of protection for areas with fragile skin tissue.

HCPC CODE E0191



MEDI-TECH
International Corporation

(800) 333-0109

www.medi-techintl.com