

MT SPANDAGE™

TUBULAR ELASTIC RETAINER NET

Not Made with Natural Rubber Latex



HCPCS Code A6457

- Replaces adhesive tape & wrap gauze
- Holds both wet & dry dressings firmly in place
- 15 sizes including MT22 with 80" circumference
- A lucite applicator is included in each box of MTF1 & MTF2 at no additional cost
- No allergic reactions. No surprises
- Eliminates skin stripping & excoriation caused by tape removal
- No "wet tape" or tape residue removal
- Equal pressure retains dressing without slippage or excessive constriction
- Dressing changes are easier, faster, & more Comfortable - less time for staff, less Stressful for patients

Product #	Description / Application
MTF1	Small finger, toe w/applicator
MTF2	Adult finger, toe w/applicator
MT01	Fingers, toes, wrists
MT02	Small hand, arm, leg, foot
MT03	Medium hand, arm, leg, foot
MT04	Large hand, arm, leg, foot
MT05	Small head, shoulder, thigh
MT06	Large head, shoulder, thigh
MT07	Small chest, back, perineum, axilla
MT08	Medium chest, back, perineum, axilla
MT09	Large chest, back, perineum, axilla
MT10	X-Large chest, back, perineum, axilla
MT11	XX-Large chest, back, perineum, axilla
MT12	XXX-Large chest, back, perineum, axilla
MT22	Custom Size (Circumference stretch to 80")

25 Yards Per Box - 10 Boxes Per Case
Also available in 10 Yards & 50 Yards Per Box.
Pre-Cuts and Pre-Cuts with Cut Outs.



(800) 333-0109

www.medi-techintl.com

MT SPANDAGE Tubular Retainer Net



Available in 15 sizes

25 yards stretched per box

10 yards stretched per box

50 yards stretched per box

Pre-Cut Sizes

Pre-Cut Sizes with Cut-Outs

SIZING OPTIONS

1. HEAD

- MT05 - X-Small
- MT06 - Small
- MT07 - Medium
- MT08 - Large



7. FOOT

- MT02 - X-Small
- MT03 - Small
- MT04 - Medium
- MT05 - Large

2. SHOULDER


- MT05 - Small
- MT06 - Medium
- MT07 - Large

4. HAND

- MT02 - Small
- MT03 - Medium
- MT04 - Large

8. THIGH

- MT05 - Small
- MT06 - Medium
- MT07 - Large
- MT08 - X-Large



3. CHEST - BACK

- AXILLA - PERINEUM**
- MT07 - Small
 - MT08 - Medium
 - MT09 - Large
 - MT10 - X-Large
 - MT11 - X-Large
 - MT12 - XX-Large
 - MT22 - Bariatric Size
 - Special Comfort Size
 - Stretches to 80 inches circumference

5. LEG

- MT02 - X-Small
- MT03 - Small
- MT04 - Medium
- MT05 - Large

9. FINGERS

- MTF1 w/applicator - Small
- MTF2 w/applicator - Medium
- MT01 - Small
- MT02 - Medium

6. TOES

- MTF1 w/applicator - Small
- MTF2 w/applicator - Medium
- MT01 - Small
- MT02 - Medium

10. ARMS

- MT02 - X-Small
- MT03 - Small
- MT04 - Medium
- MT05 - Large

Relaxed and Stretched Sizes of MT Spandage Tubular Retainer Net

Item #	Size	Relaxed	10" Max. Stretch
MTF1	Size .05	.05"	1.00"
MTF2	Size .05	.05"	1.50"
MT01	Size 1	0.625"	6.00"
MT02	Size 2	1.00"	7.50"
MT03	Size 3	1.125"	10.00"

Item #	Size	Relaxed	10" Max. Stretch
MT04	Size 4	1.25"	10.50"
MT05	Size 5	1.50"	14.50"
MT06	Size 6	2.00"	18.50"
MT07	Size 7	2.75"	27.50"
MT08	Size 8	3.00"	34.00"

Item #	Size	Relaxed	10" Max. Stretch
MT09	Size 9	3.25"	39.00"
MT10	Size 10	3.750"	40.50"
MT11	Size 11	4.50"	47.00"
MT12	Size 12	5.03"	58.00"
MT22	Size 22	7.75"	81.00"