

Made in
the USA

SPANDAGRIP™

TUBULAR ELASTIC SUPPORT BANDAGE

Not Made with Natural Rubber Latex

SpandaGrip™ stays in place without pins, clips or tape. A convenient alternative to flat wrap bandages. Easy to apply. Not bulky, which allows freedom of movement and comfort.

Delivers consistent support with a single layer.



PRODUCT INFORMATION

Perfect for the effective management of:

- Sprains and Strains
- Support following Cast Removal
- Soft-tissue Injuries
- Joints and Muscles
- General Edema
- Post-graft Pressure for Burn Patients

Provided in 11 Yards / 33 Feet lengths per box.

Order #	Description and Color	Size	Suggested Application
SAG13110	A - Natural	1-1/2"	Infant feet and arms
SAG13111	B - Natural	2-1/2"	Small hands and arms
SAG13112	C - Natural	2-3/4"	Medium arms, small ankles
SAG13113	D - Natural	3"	Large arms, med. ankles, small knees
SAG13114	E - Natural	3-1/2"	Large ankles, med. knees, small thighs
SAG13115	F - Natural	4"	Large knees, medium thighs
SAG13116	G - Natural	4-1/2"	Large thighs
SAG13117	J - Natural	6-3/4"	Small trunks
SAG13118	K - Natural	8-1/4"	Medium trunks
SAG13119	L - Natural	13"	Large trunks
SAG13120	H - Natural	5"	X-Large thighs
SAG13121	I - Natural	5-1/2"	X-Small trunks
SAG13140	B - Beige	2-1/2"	Small hands and arms
SAG13141	C - Beige	2-3/4"	Medium arms, small ankles
SAG13142	D - Beige	3"	Large arms, med. ankles, small knees
SAG13143	E - Beige	3-1/2"	Large ankles, med. knees, small thighs
SAG13144	F - Beige	4"	Large knees, medium thighs
SAG13145	G - Beige	4-1/2"	Large thighs



MEDI-TECH
International Corporation

(800) 333-0109

www.medi-techintl.com

SpandaGrip™

Tubular Elastic Support Bandage

Not Made with Natural Rubber Latex

Provides Sustained Pressure & Support for: Strains & Sprains, Soft Tissue Injuries, Support following Cast Removal, General Edema, Post-Graft Pressure for Burn Patients & Rib Cage Injuries

INSTRUCTIONS FOR USE



1. CUT A LENGTH OF SpandaGrip™ TO FIT THE LIMB OR AREA TO BE COVERED -ALLOW AN EXTRA 2-3cm (1 inch) FOR OVERLAP
2. DRAW THE BANDAGE CAREFULLY OVER THE LIMB OR AREA LIKE A STOCKING
3. A SINGLE LAYER ONLY IS RECOMMENDED
4. SIZING VARIANCE DUE TO BODY TYPE

SPANDAGRIP TUBULAR ELASTIC BANDAGE - COMPRESSION SIZING CHART

SINGLE LAYER APPLICATION

NATURAL			Low Pressure 5 - 10 mmHg		Medium Pressure 10 - 20 mmHg		High Pressure 20 - 30 mmHg	
PRODUCT NUMBER	Color Code	SIZE	Low Pressure Circumference Centimeters	Circumference Inches	Medium Pressure Circumference Centimeters	Circumference Inches	High Pressure Circumference Centimeters	Circumference Inches
SAG13110	Purple	1.50" (A)	19 - 22 cm	7.5 - 8.7 in	22 - 24 cm	8.7 - 9.4 in	24 - 27 cm	9.4 - 10.6 in
SAG13111	Olive Green	2.50" (B)	22 - 24 cm	8.7 - 9.4 in	24 - 27 cm	9.4 - 10.6 in	27 - 30 cm	10.6 - 11.8 in
SAG13112	Yellow	2.75" (C)	24 - 27 cm	9.4 - 10.6 in	27 - 30 cm	10.6 - 11.8 in	30 - 35 cm	11.8 - 13.8 in
SAG13113	Light Blue	3.00" (D)	27 - 30 cm	10.6 - 11.8 in	30 - 35 cm	11.8 - 13.8 in	35 - 44 cm	13.8 - 17.3 in
SAG13114	Brown	3.50" (E)	30 - 35 cm	11.8 - 13.8 in	35 - 44 cm	13.8 - 17.3 in	44 - 48 cm	17.3 - 18.9 in
SAG13115	Aqua	4.00" (F)	35 - 44 cm	13.8 - 17.3 in	44 - 58 cm	17.3 - 22.8 in	58 - 64 cm	22.8 - 25.2 in
SAG13116	Pink	4.50" (G)	44 - 58 cm	17.3 - 22.8 in	58 - 64 cm	22.8 - 25.2 in	64 - 72 cm	25.2 - 28.3 in
SAG13121	Dark Blue	5.50" (I)	58 - 64 cm	22.8 - 25.2 in	64 - 72 cm	25.2 - 28.3 in	72 - 86 cm	28.3 - 33.8 in
SAG13117	Orange	6.75" (J)	64 - 72 cm	25.2 - 28.3 in	72 - 86 cm	28.3 - 33.8 in	86 - 96 cm	33.8 - 37.8 in
SAG13118	Red	8.25" (K)	72 - 86 cm	28.3 - 33.8 in	86 - 96 cm	33.8 - 37.8 in	96 - 108 cm	37.8 - 42.5 in

PRODUCT NUMBER	Color Code	SIZE	Low Pressure Circumference		Medium Pressure Circumference		High Pressure Circumference	
			Centimeters	Inches	Centimeters	Inches	Centimeters	Inches
SAG13120	Black	5.00" (H)	53 - 62 cm	20.9 - 24.4 in	62 - 70 cm	24.4 - 27.6 in	70 - 83 cm	27.6 - 32.7 in
SAG13119	Light Green	13.00" (L)	105 - 130 cm	41.3 - 51.2 in	130 - 149 cm	51.2 - 58.7 in	149 - 160 cm	58.7 - 63 in

SpandaGrip™ has been designed for tissue support and provides sustained pressure. It does not provide graduated therapeutic compression. The application of SpandaGrip is solely at the discretion of the treating physician or other certified clinical personnel and should not be utilized without considering the clinical presentation of the patient and the facility's protocols and guidelines. Any information, sizing charts or other sizing instruments provided act merely as guidelines to assist in the proper application of SpandaGrip, not as recommended usage for any patient or clinical condition.



Made in
the USA

SPANDAGRIP™

TUBULAR ELASTIC SUPPORT BANDAGE

SINGLE UNITS

Not Made with Natural Rubber Latex

SpandaGrip™ stays in place without pins, clips or tape. A convenient alternative to flat wrap bandages. Easy to apply. Not bulky, which allows freedom of movement and comfort.

Delivers consistent support with a single layer.



Order #	Description and Color	Size	Suggested Application
SAG25036	B - Natural	2-1/2" x 36"	Small hands and arms
SAG23436	C - Natural	2-3/4" x 36"	Medium arms, small ankles
SAG10336	D - Natural	3" x 36"	Large arms, med. ankles, small knees
SAG35036	E - Natural	3-1/2" x 36"	Large ankles, med. knees, small thighs
SAG10436	F - Natural	4" x 36"	Large knees, medium thighs
SAG45036	G - Natural	4-1/2" x 36"	Large thighs

PRODUCT INFORMATION

Perfect for the effective management of:

- Sprains and Strains
- Support following Cast Removal
- Soft-tissue Injuries
- Joints and Muscles
- General Edema
- Post-graft Pressure for Burn Patients

36 Inches Per Unit / 12 Units Per Case



(800) 333-0109

www.medi-techintl.com

Made in
the USA

MEASURING TAPE

DESIGNED TO DETERMINE APPLICATION SIZE OF

SPANDAGRIP™
TUBULAR ELASTIC SUPPORT BANDAGE

Order #	PKG	Description	Suggested Application
SAGTAPE	12/Case	SpandaGrip™ Measuring Tape measures from 6" to 43" and/or 15cm to 109cm Color Coded by Size	Measure the largest circumference of the area to be treated with low, medium or high compression Always apply SpandaGrip™ as a SINGLE LAYER!



How to Clean the Measuring Tape:

For sanitary reasons, MTIC recommends using the measuring tape on a single patient. If your protocol indicates that it can be used on multiple patients, the tape measure can be sanitized by wiping the entire length of the tape (both sides) with an alcohol prep pad. Let it dry completely before storing it and reusing

SpandaGrip™ has been designed for tissue support and provides sustained pressure. It does not provide graduated therapeutic compression. The application of SpandaGrip™ is solely at the discretion of the treating physician or other certified clinical personnel and should not be utilized without considering the clinical presentation of the patient and the facility's protocols and guidelines. Any information, sizing charts or other sizing instruments provided act merely as guidelines to assist in the proper application of SpandaGrip™, not as recommended usage for any patient or clinical condition.



MEDI-TECH
International Corporation
(800) 333-0109

www.medi-techintl.com